


Part No. 98321: 5 ml Adapter for Use with 98316 Mounting Bracket Kit



Description:

This is inserted into a 98316 mounting bracket kit to allow it to accept 5 ml syringes. **The 98316 mounting bracket must be purchased separately.**

Includes:

Qty	Description
1	5 ml adapter for use with 98316

Required: *** Must Be Purchased Separately ***

Part No.	Description
98316	30 ml mounting bracket kit

Instructions for Use:

1. Insert the adapter into the mounting bracket kit.
2. Tighten the mounting bracket using the two black screws.